INNOSOL

Mesa 160

User's Instructions
Gebrauchsanleitung
Manuel de l'utilisateur
Gebruiksaanwijzing
Bruksanvisningar
Käyttäjän opas

Made in Finland


Innojok Oy Ltd
Helsinki

SPECTRIAL DISTRIBUTION DIAGRAMS

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Congratulations for purchasing **INNOSOL Mesa 160** bright light therapy lamp. Before using the lamp, it is good to know what effects bright light has on your body and how to use bright light therapy effectively. Many people suffer from lethargy and depression during the dark winter season. This condition is known as SAD (Seasonal Affective Disorder). Scientific research has shown that as the amount of natural light decreases, the body may physiologically respond to these "unfavourable living conditions" by adjusting certain hormone levels and biological activities. The effects of this adjustment may include:

- Intense fatigue, exhaustion and various sleep disorders
- Increased appetite leading to weight gain
- General passivity and depression
- Decreased work efficiency and loss of vitality

The strength of light on a sunny day is approximately 50 000 – 100 000 lux, whereas indoors it is as little as 100-200 lux. An effective bright light therapy session requires the eyes to be exposed to at least 2 500 lux. Via the retina and the optic nerves the light affects the brain's light sensitive region (suprachiasmatic cellkernel), which adjusts the body's physiological cycles in accordance with the rhythm of day and night. The light reduces the amount of the sleep hormone melatonin, secreted by the pineal gland. Melatonin is a key factor in the regulation of alertness and sleep cycles, but also in the prevention of SAD. Bright light therapy helps you alleviate these problems and regain vitality.

If the symptoms are severe, they should be treated under medical supervision.

THERAPY INSTRUCTIONS

Start bright light therapy before any symptoms of SAD appear, preferably soon after nights start drawing in.

Take a bright light session 3-5 times a week throughout the winter. If the symptoms are not relieved or they reappear, try taking daily sessions for 1-2 weeks. You may even begin with daily sessions or repeat this intensity whenever you feel you need it.

One session lasts for 0,5-2 hours. Longer than two hours does not intensify the effect. It is advisable to take the session between 6.00 a.m. and 6.00 p.m. but it is most effective in the morning. A session in the evening may disturb sleep. If you tend to have evening sleepiness (common to elderly people), you might want to try taking the sessions earlier in the evening in order to stay up longer. You can also divide the session into several shorter periods during the same day if it suits your schedule better.

There is no need to look straight into the light source, but you should keep it within your visual field during the whole session. You can read, write or do other things in the meantime. The effectiveness of bright light therapy correlates with the distance between your face and the light source. **Make sure that your face is no more than 105 cm away from the lamp.** You can shorten the required session time by sitting closer to the lamp.

Bright light therapy also helps speed recovery from jet-lag. A session in the evening will delay your sleep cycle, a session in the morning will have an opposite effect.

IMPORTANT NOTES ON BRIGHT LIGHT THERAPY

The lamp comes with an electronic ballast which raises the light frequency up to 30 000 Hz, where flicker is not noticeable to the eye. After the first few sessions, mild headache or eye strain may occur. This is quite normal and such side effects usually disappear in a few days. If this does not happen, stop the therapy and consult your doctor before continuing.

Not all people suffer from SAD due to lack of daylight. Exhaustion, sleep disorders and depression may have other causes.

Bright light therapy is most effective when combined with regular physical exercise and a healthy diet.

Consult your doctor before starting bright light therapy if you suffer from

- **eye illnesses (e.g. eye inflammation or hypersensitivity)**
- **mental disorders**
- **or if you use psychotropic drugs**

When you start the bright light therapy, it is useful to take notes on the effects it has on you. This helps you determine which way the therapy works best for you. The need for bright light therapy as well as its duration and frequency vary individually. Use the therapy diary on the back of this booklet.

1. Install the fluorescent tubes by following the steps 2,3,5 and 6 indicated in "Replacing the fluorescent tubes".
2. Place the lamp upright on a table using the support. The lamp is designed to not use in a bathroom, sauna or any other damp place, nor outdoors.
3. If you attach the lamp to a wall, first fasten a screw firmly in the wall and hang the lamp from the hole provided.
4. Connect the plug to a wall socket. The charge must be 220 - 240 V.
5. Use the switch in the electric cord to turn on the lamp. The light levels of the new fluorescent tubes may fluctuate during the first few hours, but this usually ceases after approximately ten hours of use.
6. Switch off the power after use.

REPLACING THE FLUORESCENT TUBES

1. When replacing the tubes, check that the new ones have similar specifications to the originals.
 2. Place the lamp on its back on a dry, clean and flat surface. Make sure that the plug is disconnected.
 3. Remove the protective shield by gently pressing on its top edge while pulling outwards. Remove the shade from the lamp.
 4. Firmly grip the base of the fluorescent tube and detach it from the lampholder by pulling it upwards. Detach the upper end of the tube from its support bracket. Remove the other tube in the same way if necessary.
 5. Install the new tube by placing the upper end in the support bracket. Insert tube's contacts into the holes in the lampholder. Firmly grip the base of the tube and push it into the lampholder until it clicks into place. Install the other tube in the same way if needed.
 6. Once the fluorescent tubes are firmly in place, check that they work before replacing the shade.
- Please do not discard tubes with normal household refuse. Dispose them appropriately so as not to endanger people or the environment.

MAINTENANCE OF THE LAMP

Keep the lamp clean and dust free. Always disconnect and allow tubes to cool down before cleaning. Use a damp cloth to clean its outer surfaces as well as the fluorescent tubes and inner surfaces. Do not use any solvents or detergents which may harm the surfaces. Do not let any water leak into the lamp.

IMPORTANT

- Do not let children play with the lamp
- If the electric cord is damaged, it must be replaced by an authorised electrician

TECHNICAL SPECIFICATIONS

- weight 4,9 kg
- length 65 cm, width 40,5 cm, depth 11 cm
- fluorescent tubes 2 x 80 W, 4 000 K
- surface luminosity 11 000 cd/m², max. 18 000 cd/m²
- electronic ballast
- charge 220 – 240 V/50 Hz
- 3,0 m electric cord with earthed plug
- opaque glare- and UV-protective shade
 - powder painted body of zinc plated steel
- adjustable leg
- wall mount option
- therapy distance 105 cm (2 500 lx)
- FI approved by the Electric Inspectorate of Finland (SGS Fimko), certified DNV0434 Class 2A

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