INNOSOL Rondo

**User's Instructions
Gebrauchsanleitung
Manuel de l'utilisateur
Gebruiksaanwijzing
Bruksanvisningar
Käyttäjän opas**

INNOJOK LTD
Helsinki, Finland
www.innojok.fi

English	2
Deutsch	4
Français	6
Nederlands.....	8
Svenska	10
Suomi	12

Congratulations for purchasing INNOSOL Rondo bright light therapy lamp. Before using the lamp, it is good to know what effects bright light has on your body and how to use bright light therapy effectively. Many people suffer from lethargy and depression during the dark winter season. This condition is known as SAD (Seasonal Affective Disorder). Scientific research has shown that as the amount of natural light decreases, the body may physiologically respond to these "unfavourable living conditions" by adjusting certain hormone levels and biological activities. The effects of this adjustment may include:

- Intense fatigue, exhaustion and various sleep disorders
- Increased appetite leading to weight gain
- General passivity and depression
- Decreased work efficiency and loss of vitality

The strength of light on a sunny day is approximately 50 000 – 100 000 lux, whereas indoors it is as little as 100-200 lux. An effective bright light therapy session requires the eyes to be exposed to at least 2 500 lux. Via the retina and the optic nerves the light affects the brain's light sensitive region (suprachiasmatic cellkernel), which adjusts the body's physiological cycles in accordance with the rhythm of day and night. The light reduces the amount of the sleep hormone melatonin, secreted by the pineal gland. Melatonin is a key factor in the regulation of alertness and sleep cycles, but also in the prevention of SAD. Bright light therapy helps you alleviate these problems and regain vitality.

If the symptoms are severe, they should be treated under medical supervision.

THERAPY INSTRUCTIONS

Start bright light therapy before any symptoms of SAD appear, preferably soon after nights start drawing in.

Take a bright light session 3-5 times a week throughout the winter. If the symptoms are not relieved or they reappear, try taking daily sessions for 1-2 weeks. You may even begin with daily sessions or repeat this intensity whenever you feel you need it.

One session lasts for 0,5-2 hours. Longer than two hours does not intensify the effect. It is advisable to take the session between 6.00 a.m. and 6.00 p.m. but it is most effective in the morning. A session in the evening may disturb sleep. If you tend to have evening sleepiness (common to elderly people), you might want to try taking the sessions earlier in the evening in order to stay up longer. You can also divide the session into several shorter periods during the same day if it suits your schedule better.

There is no need to look straight into the light source, but you should keep it within your visual field during the whole session. You can read, write or do other things in the meantime. The effectiveness of bright light therapy correlates with the distance between your face and the light source. **Make sure that your face is no more than 50 cm away from the lamp.** You can shorten the required session time by sitting closer to the lamp.

Bright light therapy also helps speed recovery from jet-lag. A session in the evening will delay your sleep cycle, a session in the morning will have an opposite effect.

IMPORTANT NOTES ON BRIGHT LIGHT THERAPY

The lamp comes with an electronic ballast which raises the light frequency up to 30 000 Hz, where flicker is not noticeable to the eye. After the first few sessions, mild headache or eye strain may occur. This is quite normal and such side effects usually disappear in a few days. If this does not happen, stop the therapy and consult your doctor before continuing. Not all people suffer from SAD due to lack of daylight. Exhaustion, sleep disorders and depression may have other causes.

Bright light therapy is most effective when combined with regular physical exercise and a healthy diet.

Consult your doctor before starting bright light therapy if you suffer from

- eye illnesses (e.g. eye inflammation or hypersensitivity)
- mental disorders
- or if you use psychotropic drugs

When you start the bright light therapy, it is useful to take notes on the effects it has on you. This helps you determine which way the therapy works best for you. The need for bright light therapy as well as its duration and frequency vary individually. Use the therapy diary on the back of this booklet.

1. Place the lamp on a table with the shade facing up. **The fluorescent tube is inside the lamp, protected in cardboard.**
2. Open the shade by turning a small coin or a screw driver in one of the notches on the outer rim of the shade.
3. Remove the cardboard around the tube and insert the tube back into the lampholder.
4. Attach the shade by pressing its rim firmly into the groove on the lamp frame.
5. If you attach the lamp to a wall, first fasten a screw firmly in the wall and hang the lamp in the hole provided.
6. Push the front edge of the metal plate (fastened in the back of the lamp) under the bent edge of the other metal plate. Place the lamp on its support leg upright on a table.
7. Connect the plug to a wall socket. The charge must be 230 V.
8. Use the switch in the cord to turn on the lamp. The light levels of the new fluorescent tubes may fluctuate during the first few hours, but this usually ceases after approximately ten hours of use.
9. Switch off the power after use.

REPLACING THE FLUORESCENT TUBE

1. When replacing the tube, check that the new one has similar specifications to the original.
2. Place the lamp on its back on a flat surface. Make sure that the plug is disconnected.
3. Detach the shade, replace the tube and attach the shade (steps 2-4 above).
4. Please do not discard tubes with normal household refuse. Dispose them appropriately so as not to endanger people or the environment.

MAINTENANCE OF THE LAMP

Keep the lamp clean and dust free. Always disconnect and allow the tube to cool down before cleaning. Use a damp cloth to clean its outer surfaces as well as the fluorescent tube and inner surfaces. Do not use any solvents or detergents which may harm the surfaces. Do not let any water leak into the lamp.

IMPORTANT

- Do not let children play with the lamp
- If the electric cord is damaged, it must be replaced by an authorised electrician

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